


TERRAWHEELS - Technical Guide



I. GENERAL CHARACTERISTICS

Product: TerraWheels is an all-terrain wheelchair for people with reduced mobility ideal for walking on soft grounds.

SUMMARY

- I. General Characteristics
- II. Product accessories
- III. Applications and installation conditions
- IV. Warranty
- V. References





Dimensions :


- Width (between armrests) :** 80 cm
- Length :** 120 cm
- Height :** 117 cm
- Weight: :** 40 kilos

Colours : Blue and white vertical stripes, yellow wheels

Materials : Its wide backrest made with a UV and seawater-resistant fabric is designed to offer a maximum comfort while sunbathing.

- Components :**
- Wheels: rubber pneumatic
 - Axle in stainless steel
 - Chassis: folding aluminium
 - Stainless trolley

- Advantages :**
- Comfortable
 - Very good stability on sand
 - **Prompt maintenance**
 - Foldable: easy to store and transport



II. PRODUCT ACCESSORIES

Back stow-away pocket:

It brings more comfort and practicality for the user.

Space to insert a parasol:

Accessrec includes a space to insert a parasol. Be careful, any parasol is included with TerraWheels.

Seatbelt:

TerraWheels is provided with a seatbelt to avoid any risk of fall when movement on sand or into water. So, users can enjoy the beach, ocean, lakes, waterparks and pools with stability and comfort.

Footrest:

TerraWheels possesses a footrest to ensure comfort on use.



III. APPLICATIONS AND INSTALLATION CONDITIONS

TerraWheels beach wheelchair is specifically dedicated to the movement of people with reduced mobility on soft grounds. Always use TerraWheels with the close assistance and supervision of an able-bodied user.

IV. WARRANTY

Accessrec Europe offers a two (2) years warranty starting at the date of your Purchase Order. This warranty only covers manufacturing defects of the wheelchair frame. Warranty does not cover fabric, armrest or pneumatic. Each guarantee limit will be assessed by Accessrec Europe.

The warranty does not cover:

- The normal use of TerraWheels
- Improper care, operation, maintenance and storage of TerraWheels
- An accident of natural disaster
- Modifications that have been made without written consent of Accessrec Europe
- Maintenance and transport costs that can be necessary if reparations are required

If you are not satisfied by the intended conditions of reparations that have been made, please contact directly Accessrec Europe. Products may not be returned without consent from Accessrec Europe. The two (2) years limited warranty is applicable to the original purchaser only.

V. REFERENCES



TerraWheels - Malibu (USA)



TerraWheels - Bordeaux, Gironde